

---

## Table Of Contents

---

Section	Description	Page #
<b>7.</b>	<b>Micro Linear Pump Module .....</b>	<b>7-2</b>
7.1	Description .....	7-2
7.2	Operation .....	7-2
7.2.1	Piston/Cylinder Set .....	7-2
7.3	Installation .....	7-3
7.4	Options .....	7-3
7.4.1	Case Material .....	7-3
7.4.2	Internal Pump Modification .....	7-3
7.4.3	Fittings .....	7-3
7.4.4	Teflon End Cap/Seal .....	7-3
7.5	Maintenance .....	7-3
7.5.1	Preventative Maintenance .....	7-3
7.5.2	O-Rings .....	7-5
7.5.3	Assembly/Disassembly Procedures .....	7-5
7.6	Problem Guide .....	7-8
7.6.1	Piston Seized In The Cylinder .....	7-8
7.7	Specifications .....	7-8
7.8	Model Number .....	7-8
7.9	Illustrated Parts Breakdown .....	7-10

## 7. MICRO LINEAR PUMP MODULE

### 7.1 DESCRIPTION (Figure 7.1)

The Micro Linear Pump Module, hereafter referred to as the Pump Module, is comprised of the following major components; a ceramic piston fabrication, a cylinder pressed into a case and intake and discharge ports designed to accept a 1/4-28 or 5/16-28 male threaded fitting. The Pump Module is within the liquid path and is designed to be detached from the Actuator Module and completely disassembled for ease of cleaning, decontamination and sterilization. The Pump Module is designed to be used in conjunction with the 40-pitch Actuator Module. Installing this Pump Module on a 20-pitch Actuator Module will damage the Pump Module.

### 7.2 OPERATION

The Pump Module is a piston/cylinder arrangement providing positive displacement. The Pump Module contains a ceramic piston that is selectively rotated and reciprocated by the Actuator Module. The piston is connected to the Actuator Module through a unilaterally flexible coupling located on the end of the Actuator Module's lead screw. The piston incorporates a flat on one end that provides valving of the pump to either the intake or discharge ports.

The piston's home position is with the piston flat aligned with the discharge port and retracted to fill the cylinder with liquid. The piston is pushed forward forcing the required amount of liquid through the discharge port. Depending on the Controller Module mode, the piston will either stop after completing a single dispense and wait until requested to dispense again or will automatically rotate to the intake port, retract to fill the chamber and then rotate back to the discharge port.

The end of the piston is never drawn back beyond the intake and discharge ports in normal operation. The piston flat allows only one port to communicate with the interior of the pump cylinder at any time. The effect of this is positive mechanical valving, eliminating the need for check valves under normal operations.

The pump, which cannot be driven by liquid pressure, essentially acts as a closed valve when the unit is not in operation.

#### 7.2.1 Piston/Cylinder Set

The piston/cylinder set is constructed of high density alumina ceramic and/or magnesium partially stabilized zirconia ceramic. The ceramics are compatible with most acids and bases. The piston/cylinder set has a clearance between the piston and cylinder wall of approximately .00005" which minimizes fluid slip.

The ceramic piston operates within the ceramic cylinder with no lubrication other than the liquid being dispensed or metered. The natural crystalline structure of the ceramic displays zero porosity ensuring zero retention and carryover of one liquid to the next.




Figure 7.1 Micro Linear Pump Module

The ceramic material's mechanical and thermal stability allows the Pump Module to be sealed by virtue of a close running clearance between the piston and the cylinder bore. This means that no compliant dynamic seals are used eliminating a part requiring frequent replacement in traditional pump designs.

### 7.3 INSTALLATION

No installation of the Pump Module is required. Refer to section 7.5.3 for assembly and disassembly of the Pump Module to the Actuator Module.

### 7.4 OPTIONS

#### 7.4.1 Case Material

The case material is Tefzel.

#### 7.4.2 Internal Pump Modification

Modifications to the piston or cylinder can be made to provide improved pumping performance for certain liquids.

#### 7.4.3 Fittings

Some of the fittings available from IVEK Corporation are listed below. Refer to Chapter 9 – Fittings and Miscellaneous Items for a complete list.

- Barb Fittings x 1/4-28  
1/16", 1/8" and 3/16" (barb size)
- Flangeless Compression Fittings  
1/4-28 (use with 1/16" or 1/8" O.D. plastic tubing)  
5/16-24 (use with 3/16" O.D. plastic tubing)

#### 7.4.4 Teflon End Cap/Seal

The Teflon end cap/seal replaces the sight glass end cap and the inner O-ring. The Teflon provides the seal for the head end of the pump cylinder. The assembly/disassembly procedures contain separate *italicized* instructions for this option and an illustrated parts breakdown provides a listing of the parts.

### 7.5 MAINTENANCE (Figure 7.2)

#### CAUTION

*Never connect or disconnect the cable from the Controller or Actuator Modules connector while power is on. Damage to the equipment may result.*

#### 7.5.1 Preventative Maintenance

The ceramic components for the pump have been designed to last for millions of repetitions without wear. Preventative maintenance includes careful handling of the piston fabrication and cylinder housing when they have been removed from the pump assembly. Always take great care when removing the piston fabrication from the cylinder and replacing the piston fabrication into the cylinder. If the cleaning procedure includes removing the Pump Module and individually cleaning separate parts, always keep the Pump Module parts together, each piston fabrication with the cylinder housing to which it was originally mated. The number on the piston fabrication should match the number on the case. Never clean in such a way that the ceramics can vibrate against each other, or chipping may result.

**CAUTION**


*Ceramic piston/cylinder sets are particularly sensitive to neglect and may seize if allowed to dry out without adequate cleaning.*

**7.5.1.1 General Applications; Routine Cleaning Procedure.**

1. Disconnect intake tubing from process liquid supply container.
2. Cycle pump in continuous mode until remaining process liquid has been purged from the Pump Module liquid path.
3. Connect the intake tubing to the cleaning liquid supply container.
4. Cycle pump in continuous mode at a high prime rate to flush the cleaning liquid through the entire liquid path.

**NOTE**

*Routine flushing with a compatible liquid after shutdown may suffice for most applications.*



**Figure 7.2 Micro Linear Pump Module Components Assembly/Disassembly**

032219 with Teflon End  
Cap/Seal, Site Glass variant not  
shown.




Figure 7.3 A-Series Linear Pump Module with Lip Seal Components Assembly/Disassembly

### 7.5.2 O-Rings (Figure 7.2)

The inner O-ring (9B) serves to seal the head end of the pumping chamber.

Over time, O-rings may lose elasticity and become deformed. Periodic Replacement of these O-rings is required. The replacement cycle is dependent on handling during assembly and disassembly in addition to the liquids being pumped.

Please contact technical support at IVEK Corporation with any questions or concerns you may have regarding the operation or maintenance of this module.

### 7.5.3 Assembly/Disassembly Procedures (Figure 7.2 & 7.3)

The Pump Module contains the following replaceable parts. Also contained in this section are the procedures for assembling and disassembling the Pump Module from the Actuator Module.


- End Cap Retainer (7)
- Sight Glass End Cap (10) or *optional Teflon End Cap/Seal (10)*
- O-Rings (9A and 9B)
- Cylinder/Case and Piston Assembly (3 and 5)

#### WARNING

*Make sure the power is OFF and all hazardous liquids have been flushed from the system prior to performing any disassembly or assembly procedures.*

#### 7.5.3.1 Pump Module (Figure 7.4)

IVEK systems are shipped with the Pump Module assembled onto the Actuator Module. The following procedures are only necessary if you received a new Pump Module or for disassembly/assembly of the Pump Module for maintenance, repair or change over.



**Figure 7.4 Micro Linear Pump Module Assembly/Disassembly**

## Disassembly:

1. Loosen (do not remove) four #10-32 socket head cap screws (3) securing Pump Module (5) to Actuator Module (1).
2. Turn Pump Module (5) slightly counterclockwise and slowly remove.

**NOTE**

*Piston (4) will stay attached to coupling (2). Do not loosen screws holding coupler to lead screw.*

3. Loosen two thumb screws (8) securing access cover (9) to Actuator Module (1).
4. Slide access cover (9) to one side.
5. Loosen (do not remove) two set screws (6) securing piston (4) in hub (2).
6. Remove piston (4).

## Assembly:

1. Loosen two thumb screws (8) securing access cover (9) to Actuator Module (1).
2. Rotate access cover (9) to one side.
3. Remove piston (4) from Pump Module (5).
4. Slide the drive pin on piston (4) into coupling (2) slot. Make sure the drive pin is fully inserted and is resting against the bottom of the slot (7).
5. Secure piston (4) in coupling (2) with two set screws (6). Tighten firmly and evenly.
6. Rotate access cover (9) on Actuator Module (1) back into place and tighten two thumb screws (8).
7. Carefully align Pump Module (5) over piston (4) and gently slide it over piston (4) and four screws (3) until it abuts Actuator Module (1). Turn slightly clockwise and secure with four #10-32 socket head cap screws (3).

**NOTE**

*Make sure the intake and discharge ports are properly orientated. (same as prior to removal)*

**7.5.3.2 End Cap Retainer, End Cap and O-rings (Figure 7.2 & 7.3 items 7, 10, 9)**

## Disassembly

1. Remove four #4-40 socket head cap screws (8) securing end cap retainer (7) to cylinder housing (5).
2. Remove end cap retainer (7).
3. Remove inner O-ring (9B), sight glass end cap (10) and outer O-ring (9A).  
*OR if optional Teflon end cap/seal purchased.*
3. *Remove Teflon end cap/seal (10) and outer O-ring (9A).*

**NOTE**

*Clean and inspect seals (9 or 10) for damage prior to assembly. (Replace if necessary)*

## Assembly

1. Position outer O-ring (9A), sight glass end cap (10) and inner O-ring (9B) into the recessed diameter of end cap retainer (7).  
*OR if optional Teflon end cap/seal purchased.*
1. *Position outer O-ring (9A) and Teflon end cap/seal (10) into the recessed diameter of end cap retainer (7).*

**NOTE**

*Install sight glass end cap with the reduced diameter facing towards the cylinder housing.*

2. Position end cap retainer (7), with installed components into the recessed diameter of cylinder/case fabrication (5).
3. Secure with four #4-40 socket head cap screws (8) and torque screws (8) to 5 in lbs.

### 7.5.3.3 Adapter Ring (Figure 7.4)

#### Disassembly

1. Remove the Pump Module from the Actuator Module as described in section 7.5.3.1 steps 1 and 2 only.
2. Remove two #8-32 socket head cap screws (6) securing cylinder/case fabrication (5) to adapter ring (4) and remove adapter ring (4).

#### NOTE

*Record orientation of the Pump Module mounting holes for reassembly.*

#### Assembly

1. Position cylinder/case fabrication (5) on adapter ring (4) and secure with two #8-32 socket head cap screws (6) and torque to 5 in lbs.
2. Assemble the Pump Module onto the Actuator Module as described in section 7.5.3.1 step 7.

### 7.5.3.4 Cylinder/Case and Piston (Figure 7.2 Items 3 and 5)

#### NOTE

*The piston and cylinder are a matched set. If either the cylinder/case or piston needs replacing, both parts must be replaced.*

#### Disassembly

1. Remove Pump Module (5) from Actuator Module (1) as described in section 7.5.3.1.
2. Remove two #8-32 socket head cap screws (6) securing cylinder/case fabrication (5) to adapter ring (4) and remove adapter ring (4).
3. Remove four screws (8) securing end cap retainer (7) to cylinder housing (5).
4. Remove end cap retainer (7).
5. Remove inner O-ring (9B), sight glass end cap (10) and outer O-ring (9A).  
*OR if optional Teflon end cap/seal purchased.*
5. *Remove Teflon end cap/seal (10) and outer O-ring (9A).*

#### NOTE

*Clean and inspect seals (9 and 10) prior to assembly. (Replace if necessary)*

#### Assembly

1. Position outer O-ring (9A), sight glass end cap (10) and inner O-ring (9B) into the recessed diameter of end cap retainer (7).  
*OR if optional Teflon end cap/seal purchased.*
1. *Position outer O-ring (9A) and Teflon end cap/seal (10) into the recessed diameter of end cap retainer (7).*

#### NOTE

*Install sight glass end cap with the reduced diameter facing towards the cylinder housing.*



2. Position end cap retainer (7), with installed components into the recessed diameter of cylinder/case fabrication (5).
3. Secure with four #4-40 socket head cap screws (8) and torque screws (8) to 5 in lbs.
4. Assemble the Pump Module onto the Actuator Module as described in section 7.5.3.1.

## 7.6 PROBLEM GUIDE

### 7.6.1 Piston Seized In The Cylinder (Figure 7.4)

If the piston seizes in the cylinder perform the following steps.

#### **CAUTION**

*DO NOT TRY TO FORCE THE PISTON FREE!*

*Damage to the piston/cylinder set or Actuator Module may occur.*

1. Loosen two thumb screws (8) securing access cover (9) to Actuator Module (1).
2. Slide access cover (9) to one side.
3. Remove four bolts (3) securing Pump Module (5) to Actuator Module (1).
4. Turn the Pump Module to gain access to the two screws (6) securing piston (4) in hub (2) and loosen (do not remove) two screws (6).
5. Carefully remove Pump Module (5) from Actuator Module (1).
6. Soak the whole assembly in a liquid compatible with the materials and process liquids.
7. After soaking, try removing the piston from the cylinder by applying a light torque to the piston using only your fingers (no tools).

If the aforementioned procedures fail, contact IVEK for technical help. It may be necessary to ship the Pump Module back to the factory. Provide a note describing, in detail, what conditions caused the seizure. It may also be necessary to return the Controller and Actuator Modules along with the Pump Module should realignment of the Pump Module drive components be required.

Table 7.1 contains a list of possible problems, causes and solutions for the Pump Module.

## 7.7 SPECIFICATIONS

Table 7.2 lists the volumetric output of the different size Pump Modules and different Controller Modules. Refer to the Title Page section of this manual for the Pump Module size provided with your system. Use the **Resolution mstep (µl)** column to determine the resolution for Microstep Controller Modules and the **Resolution (µl)** column for all other Controller Modules.

**Table 7.2 Volumetric Output Of Micro Linear Pump Modules**

Size	Chamber Capacity (µl)	Resolution (µl)	Resolution mstep (µl)
4A	20	0.010	0.002
3A	50	0.025	0.005
2A	100	0.050	0.010
1A	200	0.100	0.020

## 7.8 MODEL NUMBER

The model number provides important information about the specifics of your Pump Module. Refer to this number when calling IVEK Technical support. The model number for your Pump Module is located in the Title Page section of this manual.

<b>Sight Glass End Cap</b>	<b>032127 -</b>	<b>#</b>	<b>#</b>	<b>##</b>
<b>Teflon End Cap</b>	<b>032145 -</b>	<b>#</b>	<b>#</b>	<b>##</b>

**End Cap Retainer Material**

- 1 - White Delrin
- 2 - 303 Stainless Steel

**O-Ring Material**

- 1 - Buna-N
- 2 - Ethylene Propylene
- 3 - Kalrez
- 4 - Polyurethane
- 5 - Silicone
- 6 - Teflon
- 7 - Teflon Encapsulated Silicone
- 8 - Viton
- 9 - Teflon Encapsulated Viton

**Pump Case / Ceramics**

- 11 - 1/4 - 28, 4A HIP-YTZP/Mag Zirc
- 21 - 5/16 - 24, 4A HIP-YTZP/Mag Zirc
- 1A - 1/4 - 28, 4A HIP-YTZP/Mag Zirc, W/Gland
- 2A - 5/16 - 24, 4A HIP-YTZP/Mag Zirc, W/Gland
- 12 - 1/4 - 28, 4A Mag Zirc
- 22 - 5/16 - 24, 4A Mag Zirc
- 1B - 1/4 - 28, 4A Mag Zirc, W/Gland
- 2B - 5/16 - 24, 4A Mag Zirc, W/Gland
- 13 - 1/4 - 28, 3A HIP-YTZP/Mag Zirc
- 23 - 5/16 - 24, 3A HIP-YTZP/Mag Zirc
- 1C - 1/4 - 28, 3A HIP-YTZP/Mag Zirc, W/ Gland
- 2C - 5/16 - 24, 3A HIP-YTZP/Mag Zirc, W/Gland
- 31 - 1/4 - 28, 3A HIP-YTZP/Alumina
- 42 - 5/16 - 24, 3A HIP-YTZP/Alumina
- 3A - 1/4 - 28, 3A HIP-YTZP/Alumina, W/Gland
- 4A - 5/16 - 24, 3A HIP-YTZP/Alumina, W/Gland
- 14 - 1/4 - 28, 3A Mag Zirc
- 24 - 5/16 - 24, 3A Mag Zirc
- 1D - 1/4 - 28, 3A Mag Zirc, W/Gland
- 2D - 5/16 - 24, 3A Mag Zirc, W/Gland
- 15 - 1/4 - 28, 2A Mag Zirc
- 25 - 5/16 - 24, 2A Mag Zirc
- 1E - 1/4 - 28, 2A Mag Zirc, W/Gland
- 2E - 5/16 - 24, 2A Mag Zirc, W/Gland
- 16 - 1/4 - 28, 2A Alumina
- 26 - 5/16 - 24, 2A Alumina
- 1F - 1/4 - 28, 2A Alumina, W/Gland
- 2F - 5/16 - 24, 2A Alumina, W/Gland
- 17 - 1/4 - 28, 1A Mag Zirc
- 27 - 5/16 - 24, 1A Mag Zirc
- 1G - 1/4 - 28, 1A Mag Zirc, W/Gland
- 2G - 5/16 - 24, 1A Mag Zirc, W/Gland
- 18 - 1/4 - 28, 1A Alumina
- 28 - 5/16-28, 1A Alumina
- 1H - 1/4-28, 1A Alumina, W/Gland
- 2H - 5/16-28, 1A Alumina, W/Gland
- 19 - 1/4-28, 1A HEX SA SIC
- 29 - 5/16-24, 1A HEX SA SIC
- 1J - 1/4-28, 1A HEX SA SIC, W/Gland
- 2J - 5/16-24, 1A HEX SA SIC, W/Gland

032219 -

# # # ## # # #

**Pump Case**

- 1 - Pump Case, A-Series, For Lip Seal, White Delrin, 1/4-28 Ports
- 2 - Pump Case, A-Series, For Lip Seal, Tefzel, 1/4-28 Ports
- 3 - Pump Case, A-Series, For Lip Seal, White Delrin, 5/16-24 Ports
- 4 - Pump Case, A-Series, For Lip Seal, Tefzel, 5/16-24 Ports
- 5 - Pump Case, A-Series, For Lip Seal, Stainless Stl, 1/4-28 Ports
- 6 - Pump Case, A-Series, For Lip Seal, Stainless Stl, 5/16-24 Ports

**End Cap Retainer**

- 1 - Retainer, End Cap; White Delrin
- 2 - Retainer, End Cap; Tefzel
- 3 - Retainer, End Cap; Stainless Steel

**Base, Pump Case**

- 1 - Base, Micro Pump Case, Lip Seal; 4A
- 2 - Base, Micro Pump Case, Lip Seal; 3A
- 3 - Base, Micro Pump Case, Lip Seal; 2A
- 4 - Base, Micro Pump Case, Lip Seal; 1A

**Ceramic Set**

- 41 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 4A, Hip Ytzc/Magzirc
- 42 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 4A, Mag Zirc/Magzirc
- 31 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 3A, Hip Ytzc/Magzirc
- 32 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 3A, Hip Ytzc/Alumina
- 33 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 3A, Mag Zirc/Mag Zirc
- 34 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 3A, Mag Zirc/Alumina
- 21 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 2A, Mag Zirc/Mag Zirc
- 22 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 2A, Mag Zirc/Alumina
- 23 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 2A, Alumina/Mag Zirc
- 24 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 2A, Alumina/ Alumina
- 11 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 1A, Mag Zirc/Mag Zirc
- 12 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 1A, Mag Zirc/Alumina
- 13 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 1A, Alumina/Mag Zirc
- 14 - Piston/Cylinder Set, A-Series, W/Lip Seal C'bore; 1A, Alumina/Alumina

**Lip Seal**

- 1 - Lip Seal, Spring Energized, PTFE/Stainless Stl Spring; 4A
- 2 - Lip Seal, Spring Energized, UHMW/ Elgiloy Spring; 4A
- 3 - Lip Seal, Spring Energized, PTFE/Stainless Stl Spring; 3A
- 4 - Lip Seal, Spring Energized, UHMW/ Elgiloy Spring; 3A
- 5 - Lip Seal, Spring Energized, PTFE/Stainless Stl Spring; 2A
- 6 - Lip Seal, Spring Energized, UHMW/ Elgiloy Spring; 2A
- 7 - Lip Seal, Spring Energized, PTFE/Stainless Stl Spring; 1A
- 8 - Lip Seal, Spring Energized, UHMW/ Elgiloy Spring; 1A

032219 - # # # ## # # #

**O-Ring Material**

- 11 - O-Ring, -011, Buna-N
- 12 - O-Ring, -011, Buna-N
- 21 - O-Ring, -011, Ethylene Propylene
- 22 - O-Ring, -011, Ethylene Propylene
- 31 - O-Ring, -011, Kalrez
- 32 - O-Ring, -011, Kalrez
- 41 - O-Ring, -011, Polyurethane
- 42 - O-Ring, -011, Polyurethane
- 51 - O-Ring, -011, Silicone
- 52 - O-Ring, -011, Silicone
- 61 - O-Ring, -011, Teflon
- 62 - O-Ring, -011, Teflon
- 71 - O-Ring, -011, Teflon Encapsulated Silicone
- 72 - O-Ring, -011, Teflon Encapsulated Silicone
- 81 - O-Ring, -011, Viton
- 82 - O-Ring, -011, Viton
- 91 - O-Ring, -011, Teflon Encapsulated Viton
- 92 - O-Ring, -011, Teflon Encapsulated Viton

**Cylinder End Cap**

- 1 - End Cap, Cylinder, Teflon
- 2 - End Cap, Cylinder, Sight Glass

7.9 ILLUSTRATED PARTS BREAKDOWN

The illustrated parts breakdown (Figure 7.5) contains replacement parts for the Micro Linear Pump Module.

Table 7.1 Common Operational Problems And Solutions		
PROBLEM	PROBABLE CAUSE	POSSIBLE SOLUTION
Air evident in discharge line.	Loose/Damaged Fitting	Tighten/replace fittings.
	Loose/Damaged End Cap Seals	Clean, inspect and replace if necessary.
	Cavitation	Increase inlet tubing size or reduce Pump Module speed.
Piston seizing	Particulate materials entrapped between piston and cylinder.	Disassemble Pump Module and clean all wetted surfaces.
Fluid leaks	Improper seated or worn end cap retainer, sight glass, end cap or O-ring.	Disassemble Pump Module and clean all wetted surfaces, inspect components and replace if necessary.
	Loose/damaged fitting or tubing.	Inspect and replace if necessary.
Fluid not moving in tubing when priming, dispensing or metering.	Pump Module not properly oriented on Actuator Module.	Make certain that port holes in Pump Module line up with flat on piston.
Fluid syphoning when Pump Module is inactive.	Loose Seals	Remove, inspect for scratches or damage, replace if necessary
	High input/output pressure differential	Reduce input or output pressure
	Coupling position has moved causing port communication.	Remove end cap, actuate piston forward, measure the space between the end of the piston and the end of the cylinder. If greater than 0.020", call IVEK Technical Support
	Tubing Adaptor fitting(s) are loose	Tighten fittings





Figure 7.5 Micro Linear Pump Module (Sheet 1 of 9)



	PART NUMBER	DESCRIPTION	UNITS PER ASSY
	032127-####	Micro Linear Pump Module/Sight Glass End Cap	1
<b>Model Dwg</b>			
<b>#</b>	<b>Index #</b>	<b>Part #</b>	<b>Description</b>
<b>Tab</b>			<b>Qty</b>
	2	022044	End Cap, Sight Glass
	1	182011-C04088	Screw, Socket Head Cap, #4-40 x .88" Long
	5	182011-C08088	Screw, Socket Head Cap, #8-32 x .88" Long
	8	001576	Adaptor Ring (optional)
032127-	#	##	<b>END CAP RETAINER MATERIAL</b>
1	6	102171-001	End Cap Retainer; White Delrin
2	6	102171-003	End Cap Retainer; 303 Stainless Steel
032127-	#	##	<b>O-RING MATERIAL</b>
1	3	142294-01101	O-Ring ; -11 Buna-N
2	3	142294-01102	" ; -11 Ethylene Propylene
3	3	142294-01103	" ; -11 Kalrez
4	3	142294-01104	" ; -11 Polyurethane
5	3	142294-01105	" ; -11 Silicone
6	3	142294-01106	" ; -11 Teflon
7	3	142294-01107	" ; -11 Teflon Encapsulated Silicone
8	3	142294-01108	" ; -11 Viton
9	3	142294-01109	" ; -11 Teflon Encapsulated Viton
2127-	#	###	<b>PUMP CASE/ CERAMICS</b>
11	7	032124-110	Micro Linear Pump/Case Fab; 1/4-28, 4A HIP-YTZP/Mag Zirc
21	7	032124-210	" ; 5/16-24, 4A HIP-YTZP/Mag Zirc
1A	7	032124-3A0	" ; 1/4-28, 4A HIP-YTZP/Mag Zirc, W/Gland
2A	7	032124-4A0	" ; 5/16-24, 4A HIP-YTZP/Mag Zirc, W/Gland
12	7	032124-120	" ; 1/4-28, 4A Mag Zirc
22	7	032124-220	" ; 5/16-24, 4A Mag Zirc
1B	7	032124-3B0	" ; 1/4-28, 4A Mag Zirc, W/Gland
2B	7	032124-4B0	" ; 5/16-24, 4A Mag Zirc, W/Gland
13	7	032124-130	" ; 1/4-28, 3A HIP-YTZP/Mag Zirc
23	7	032124-230	" ; 5/16-24, 3A HIP-YTZP/Mag Zirc
1C	7	032124-3C0	" ; 1/4-28, 3A HIP-YTZP/Mag Zirc W/Gland
2C	7	032124-4C0	" ; 5/16-24, 3A HIP-YTZP/Mag Zirc W/Gland
31	7	032124-111	" ; 1/4-28, 3A HIP-YTZP/Alumina
42	7	032124-211	" ; 5/16-24, 3A HIP-YTZP/Alumina
3A	7	032124-3A1	" ; 1/4-28, 3A HIP-YTZP/Alumina, W/Gland
4A	7	032124-4A1	" ; 5/16-24, 3A HIP-YTZP/Alumina, W/Gland
14	7	032124-140	" ; 1/4-28, 3A Mag Zirc
24	7	032124-240	" ; 5/16-24, 3A Mag Zirc
1D	7	032124-3D0	" ; 1/4-28, 3A Mag Zirc, W/Gland
2D	7	032124-4D0	" ; 5/16-24, 3A Mag Zirc, W/Gland

Figure 7.4 Micro Linear Pump Module (Sheet 3 of 9)



	PART NUMBER	DESCRIPTION	UNITS PER ASSY	
	032127-####	Micro Linear Pump Module/Sight Glass End Cap	1	
<b>Model Dwg</b>				
<b>#</b>	<b>Index</b>	<b>Part</b>		
<b>Tab</b>	<b>#</b>	<b>#</b>	<b>Description</b>	
			<b>Qty</b>	
032127-	#	###	<b>PUMP CASE/ CERAMICS</b>	
15	7	032124-150	Micro Linear Pump/Case Fab; 1/4-28, 2A Mag Zirc	1
25	7	032124-250	" ; 5/16-24, 2A Mag Zirc	1
1E	7	032124-3E0	" ; 1/4-28, 2A Mag Zirc W/Gland	1
2E	7	032124-4E0	" ; 5/16-24, 2A Mag Zirc W/Gland	1
16	7	032124-160	" ; 1/4-28, 2A Alumina	1
26	7	032124-260	" ; 5/16-24, 2A Alumina	1
1F	7	032124-3F0	" ; 1/4-28, 2A Alumina, W/Gland	1
2F	7	032124-4F0	" ; 5/16-24, 2A Alumina, W/Gland	1
17	7	032124-170	" ; 1/4-28, 1A Mag Zirc	1
27	7	032124-270	" ; 5/16-24, 1A Mag Zirc	1
1G	7	032124-3G0	" ; 1/4-28, 1A Mag Zirc, W/Gland	1
2G	7	032124-4G0	" ; 5/16-24, 1A Mag Zirc, W/Gland	1
18	7	032124-180	" ; 1/4-28, 1A Alumina	1
28	7	032124-280	" ; 5/16-24, 1A Alumina	1
1H	7	032124-3H0	" ; 1/4-28, 1A Alumina, W/Gland	1
2H	7	032124-4H0	" ; 5/16-24, 1A Alumina, W/Gland	1
19	7	032124-190	" ; 1/4-28, 1A HEX SA SIC	1
29	7	032124-290	" ; 5/16-24, 1A HEX SA SIC	1
1J	7	032124-3J0	" ; 1/4-28, 1A HEX SA SIC, W/Gland	1
2J	7	032124-4J0	" ; 5/16-24, 1A HEX SA SIC, W/Gland	1

Figure 7.4 Micro Linear Pump Module (Sheet 4 of 9)

		PART NUMBER	DESCRIPTION	UNITS PER ASSY
		032145-####	Micro Linear Pump Module/Teflon End Cap	1
<b>Model Dwg</b>				
#	Index	Part	Description	Qty
Tab	#	#		
	2	102104	Teflon End Cap	1
	1	182011-C04088	Screw, Socket Head Cap, #4-40 x .88" Long	4
	3	182011-C08088	Screw, Socket Head Cap, #8-32 x .88" Long	2
032127- # # ## <b>END CAP RETAINER MATERIAL</b>				
	1	6 102171-001	End Cap Retainer; White Delrin	1
	2	6 102171-003	End Cap Retainer; 303 Stainless Steel	1
032145- # # ## <b>O-RING MATERIAL</b>				
	1	3 142294-01101	O-Ring; -11 Buna-N	1
	2	3 142294-01102	" ; -11 Ethylene Propylene	1
	3	3 142294-01103	" ; -11 Kalrez	1
	4	3 142294-01104	" ; -11 Polyurethane	1
	5	3 142294-01105	" ; -11 Silicone	1
	7	3 142294-01107	" ; -11 Teflon Encapsulated Silicone	1
	8	3 142294-01108	" ; -11 Viton	1
	9	3 142294-01109	" ; -11 Teflon Encapsulated Viton	1
032145- # # ## <b>PUMP CASE/ CERAMICS</b>				
	11	7 032124-110	Micro Linear Pump/Case Fab; 1/4-28, 4A HIP-YTZP/Mag Zirc	1
	21	7 032124-210	" ; 5/16-24, 4A HIP-YTZP/Mag Zirc	1
	1A	7 032124-3A0	" ; 1/4-28, 4A HIP-YTZP/Mag Zirc, W/Gland	1
	2A	7 032124-4A0	" ; 5/16-24, 4A HIP-YTZP/Mag Zirc, W/Gland	1
	12	7 032124-120	" ; 1/4-28, 4A Mag Zirc	1
	22	7 032124-220	" ; 5/16-24, 4A Mag Zirc	1
	1B	7 032124-3B0	" ; 1/4-28, 4A Mag Zirc, W/Gland	1
	2B	7 032124-4B0	" ; 5/16-24, 4A Mag Zirc, W/Gland	1
	13	7 032124-130	" ; 1/4-28, 3A HIP-YTZP/Mag Zirc	1
	23	7 032124-230	" ; 5/16-24, 3A HIP-YTZP/Mag Zirc	1
	1C	7 032124-3C0	" ; 1/4-28, 3A HIP-YTZP/Mag Zirc W/Gland	1
	2C	7 032124-4C0	" ; 5/16-24, 3A HIP-YTZP/Mag Zirc W/Gland	1
	31	7 032124-111	" ; 1/4-28, 3A HIP-YTZP/Alumina	1
	42	7 032124-211	" ; 5/16-24, 3A HIP-YTZP/Alumina	1
	3A	7 032124-3A1	" ; 1/4-28, 3A HIP-YTZP/Alumina, W/Gland	1
	4A	7 032124-4A1	" ; 5/16-24, 3A HIP-YTZP/Alumina, W/Gland	1
	14	7 032124-140	" ; 1/4-28, 3A Mag Zirc	1
	24	7 032124-240	" ; 5/16-24, 3A Mag Zirc	1
	1D	7 032124-3D0	" ; 1/4-28, 3A Mag Zirc, W/Gland	1
	2D	7 032124-4D0	" ; 5/16-24, 3A Mag Zirc, W/Gland	1

Figure 7.4 Micro Linear Pump Module (Sheet 5 of 9)

	PART NUMBER	DESCRIPTION	UNITS PER ASSY
	032145-####	Micro Linear Pump Module/Teflon End Cap	1
<b>Model #</b>	<b>Dwg Index #</b>	<b>Part #</b>	<b>Description</b>
<b>Tab</b>	<b>#</b>	<b>#</b>	<b>Qty</b>
032145-	# #	###	<b>PUMP CASE/ CERAMICS</b>
15	7	032124-150	Micro Linear Pump/Case Fab; 1/4-28, 2A Mag Zirc
16	7	032124-160	" ; 1/4-28, 2A Alumina
25	7	032124-250	" ; 5/16-24, 2A Mag Zirc
1E	7	032124-3E0	" ; 1/4-28, 2A Mag Zirc W/Gland
2E	7	032124-4E0	" ; 5/16-24, 2A Mag Zirc W/Gland
26	7	032124-260	" ; 5/16-24, 2A Alumina
1F	7	032124-3F0	" ; 1/4-28, 2A Alumina, W/Gland
2F	7	032124-4F0	" ; 5/16-24, 2A Alumina, W/Gland
17	7	032124-170	" ; 1/4-28, 1A Mag Zirc
27	7	032124-270	" ; 5/16-24, 1A Mag Zirc
1G	7	032124-3G0	" ; 1/4-28, 1A Mag Zirc, W/Gland
2G	7	032124-4G0	" ; 5/16-24, 1A Mag Zirc, W/Gland
18	7	032124-180	" ; 1/4-28, 1A Alumina
28	7	032124-280	" ; 5/16-24, 1A Alumina
1H	7	032124-3H0	" ; 1/4-28, 1A Alumina, W/Gland
2H	7	032124-4H0	" ; 5/16-24, 1A Alumina, W/Gland
19	7	032124-190	" ; 1/4-28, 1A HEX SA SIC
29	7	032124-290	" ; 5/16-24, 1A HEX SA SIC
1J	7	032124-3J0	" ; 1/4-28, 1A HEX SA SIC, W/Gland
2J	7	032124-4J0	" ; 5/16-24, 1A HEX SA SIC, W/Gland

Figure 7.4 Micro Linear Pump Module (Sheet 6 of 9)

	PART NUMBER	DESCRIPTION	UNITS PER ASSY
	032219-#####	Pump Module, A-Series, Linear, with Lip Seal	1
<b>Model Dwg</b>			
<b>#</b>	<b>Index</b>	<b>Part</b>	<b>Qty</b>
<b>Tab</b>	<b>#</b>	<b>Description</b>	
	1	182011-C04100	4
	5	182011-C08088	2
		Screw, Socket Head Cap, 18-8 Inch; #4-40 x 1.0" Long	
		Screw, Socket Head Cap, 18-8 Inch; #8-32 x .88" Long	
032219-	# # # # # #	<b>LABEL</b>	
1			1
3	12	242116	1
		Label, Cylinder Identification, Pump Module, Style B	
032219-	# # # # # #	<b>PUMP CASE</b>	
1	2	072141-001	1
2	2	072141-002	1
3	2	072141-003	1
4	2	072141-004	1
5	2	102264-001	1
6	2	102264-002	1
		Pump Case, A-Ser., for Lip Seal, Wht. Delrin, 1/4-28 Ports	
		Pump Case, A-Ser., for Lip Seal, Tefzel, 1/4-28 Ports	
		Pump Case, A-Ser., for Lip Seal, Wht Delrin, 5/16-24 Ports	
		Pump Case, A-Ser., for Lip Seal, Tefzel, 5/16-24 Ports	
		Pump Case, A-Ser., for Lip Seal, SST, 1/4-28 Ports	
		Pump Case, A-Ser., for Lip Seal, SST, 5/16-24 Ports	
032219-	# # # # # #	<b>END CAP RETAINER MATERIAL</b>	
1	4	102171-001	1
2	4	102171-002	1
3	4	102171-003	1
		Retainer, End Cap; White Delrin	
		Retainer, End Cap; Tefzel	
		Retainer, End Cap; Stainless Steel	
032219-	# # # # # #	<b>BASE/PUMP CASE/LIP SEAL/SIZE</b>	
1		032147-001	
2	3	032147-002	
3		032147-003	
4		032147-004	
		Base, Micro Pump Case, Lip Seal; 4A	
		Base, Micro Pump Case, Lip Seal; 3A	
		Base, Micro Pump Case, Lip Seal; 2A	
		Base, Micro Pump Case, Lip Seal; 1A	

Figure 7.4 Micro Linear Pump Module (Sheet 7 of 9)

	PART NUMBER	DESCRIPTION	UNITS PER ASSY
	032219-#####	Pump Module, A-Series, Linear, with Lip Seal	1
<b>Model Dwg</b>			
<b>#</b>	<b>Index</b>	<b>Part</b>	<b>Qty</b>
<b>Tab</b>	<b>#</b>	<b>#</b>	
<b>PISTON/CYLINDER SETS</b>			
032219-	# # #	### # # #	
41	7	022383-21442	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 4A, HIP-YTZP/Mag Zirc 1
42	7	022383-22442	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 4A, Mag Zirc/Mag Zirc 1
31	7	022383-21332	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 3A, HIP-YTZP /Mag Zirc 1
32	7	022383-21333	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 3A, HIP-YTZP /Alumina 1
33	7	022383-22332	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 3A, Mag Zirc/Mag Zirc 1
34	7	022383-22333	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 3A, Mag Zirc/Alumina 1
21	7	022383-22222	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 2A, Mag Zirc/Mag Zirc 1
22	7	022383-22223	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 2A, Mag Zirc/Alumina 1
23	7	022383-23222	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 2A, Alumina/Mag Zirc 1
24	7	022383-23223	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 2A, Alumina/ Alumina 1
11	7	022383-22112	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 1A, Mag Zirc/Mag Zirc 1
12	7	022383-22113	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 1A, Mag Zirc/Alumina 1
13	7	022383-23112	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 1A, Alumina/Mag Zirc 1
14	7	022383-23113	Piston/Cylinder Set, A-Ser., w/Lip Seal C'bore; 1A, Alumina/Alumina 1
<b>LIP SEAL MATERIAL</b>			
032219-	# # #	### # #	
1	8	142378-005	Lip Seal, Spring Energized, PTFE/SST Spring; 4A 1
2	8	142669-005	" , UHMW/ Elgiloy Spring; 4A 1
3	8	142378-006	" , PTFE/Stainless Stl Spring; 3A 1
4	8	142669-006	" , UHMW/ Elgiloy Spring; 3A 1
5	8	142378-001	" , PTFE/Stainless Stl Spring; 2A 1
6	8	142669-001	" , UHMW/ Elgiloy Spring; 2A 1
7	8	142378-007	" , PTFE/Stainless Stl Spring; 1A 1
8	8	142669-007	" , UHMW/ Elgiloy Spring; 1A 1
<b>END CAP RETAINER MATERIAL</b>			
032219-	# # #	### # #	
1	6	102104	End Cap, Cylinder, Teflon 1
2	6	022044	End Cap, Cylinder, Sight Glass 1

Figure 7.4 Micro Linear Pump Module (Sheet 8 of 9)

		PART NUMBER	DESCRIPTION	UNITS PER ASSY
		032219-#####	Pump Module, A-Series, Linear, with Lip Seal	1
<b>Model Dwg</b>				
<b>#</b>	<b>Index</b>	<b>Part</b>	<b>Description</b>	<b>Qty</b>
<b>Tab</b>	<b>#</b>	<b>#</b>	<b>#</b>	
032219-	#	#	#	#
11	9	142294-01101	O-Ring; -11 Buna-N	1
12	9	142294-01101	" ; -11 Buna-N	2
21	9	142294-01102	" ; -11 Ethylene Propylene	1
22	9	142294-01102	" ; -11 Ethylene Propylene	2
31	9	142294-01103	" ; -11 Kalrez	1
32	9	142294-01103	" ; -11 Kalrez	2
41	9	142294-01104	" ; -11 Polyurethane	1
42	9	142294-01104	" ; -11 Polyurethane	2
51	9	142294-01105	" ; -11 Silicone	1
52	9	142294-01105	" ; -11 Silicone	2
61	9	142294-01106	" ; -11 Teflon	1
62	9	142294-01106	" ; -11 Teflon	2
71	9	142294-01107	" ; -11 Teflon Encapsulated Silicone	1
72	9	142294-01107	" ; -11 Teflon Encapsulated Silicone	2
81	9	142294-01108	" ; -11 Viton	1
82	9	142294-01108	" ; -11 Viton	2
91	9	142294-01109	" ; -11 Teflon Encapsulated Viton	1
92	9	142294-01109	" ; -11 Teflon Encapsulated Viton	2

Figure 7.4 Micro Linear Pump Module (Sheet 9 of 9)