

DIGISPENSE 1000* CONTROLLED SUBSTANCE DISPENSING SYSTEM

FEATURES AND BENEFITS

- High-speed brushless DC motor & controller
- Whisper-quiet motor
- Featuring Positive Displacement Ceramic Pumps
- Easy to control from the USB or RS232 interface
- Factory-made tubing kits with bottle cap connections
- Precision volumetric flask for accuracy checks
- Factory pre-calibrated motor/base pump module
- NEW all-in-one Controlled Substance Dispensing Station option: portable, secure & easy clean
- CE certified

The IVEK Digispense 1000* is a precise, high speed, low volume, positive displacement liquid dispensing system. The DigiSpense 1000 is capable of extremely high accuracy and repeatability and is often used for dosing pharmaceuticals, where inventory control and accuracy are critical. The wetted components include 99.8% pure Alumina ceramic piston/cylinder modules. These components are mated together to a clearance of 120 millionths of an inch. With only one moving part in the fluid path, systems operate for many years without measurable wear. A USB and RS232 interface port on the controller offers our customers a turnkey solution supported by most addiction treatment software.

IVEK Corporation has been manufacturing precision positive displacement dispensing systems specifically for use in addiction management programs for over 30 years. These systems are noted for their reliability, low maintenance and high precision fluid handling and are placed in treatment clinics worldwide.

*Enhanced replacement for IVEK's DIGISPENSE 10/RS232



NEW Digispense 1000 Controller
*enhanced replacement for the Digispense 10/RS232



single-end motor base module
7cm x 12.8cm x 10cm



NEW all-in-one portable, secure, easy clean,
controlled substance dispensing station option available



IVEK CORPORATION
10 FAIRBANKS ROAD NORTH SPRINGFIELD VERMONT USA 05150
TEL (01) 802-886-2238 FAX (01) 802-886-8274
TOLL FREE IN NORTH AMERICA 800-356-4746

